



Strapless N-95 Face Respirator

**The Protection You Need,
The Comfort You Deserve**



No more sore ears.
No more bruised faces.
No more foggy glasses and face shields.



Innovative Mask Technology

The multi-layered mask filters out 95% of 0.3-micron particles in the air while allowing the wearer to breathe freely. The streamlined design and flat profile make less bulky and easier to store than traditional masks.

The NIOSH-certified strapless N95 face mask is unlike any mask you've worn. Unlike traditional N95 masks, this strapless mask creates an **airtight seal** using a medical-grade adhesive, eliminating painful issues caused by tight elastic straps and ties and metal nose clips.

Country of Origin: USA Manufactured in Ohio



Interested in learning more? Visit our website: www.alliantbiotech.com



User Benefits

The outer layer provides strength and form to the mask.

The N95 filter layer provides electrostatic particle and germ filtration down to 0.3 microns.

The inner layer provides gentle comfort against your face.

An adhesive layer is a medical-grade adhesive located around the mask's perimeter, contouring to the face to provide an airtight seal.



Half the Cost and Ready to Ship

The self-adhering N95 masks are **half the cost** of traditional N95 masks, and we have masks in our warehouse, **ready to ship!** We will send samples to qualified companies and set up a meeting to show you how to use the mask properly.

Contact us to order yours today!

Part Number	Product Description	NIOSH Model Number
N1906S	Small Face Mask ACY (7in x 5in)	N1901S
N1906L	Large Face Mask ACY (8in x 6in)	N1901L
N1906XL	Xlarge Face Mask ACY (9in x 7in)	N1901XL

Global Safety First, LLC. NIOSH Approval TC#84A-8133
*See full IFU for proper application

Interested in learning more? Visit our website: www.alliantbiotech.com