



MRI (Magnetic Resonance Imaging) Safety Declaration

September 1, 2022

This is a Safety Declaration for anyone receiving a Magnetic Resonance Imaging exam (MRI). There is currently a warning informing anyone. The U.S. Food and Drug Administration (FDA) is informing health care providers and patients that they may be injured if they wear face masks (surgical masks, nonsurgical masks, and respirators) with metal parts and/or metal coatings during a Magnetic Resonance Imaging (MRI) exam. The metal parts attached to many masks, including nose pieces or nose clips, wires, nanoparticles (ultrafine particles), or antimicrobial coating that may contain metal (such as silver or copper), can become hot, burning the patient during an MRI.

However, the FDA is recommending patients wear face masks with no metal during MRIs. Alliant Biotech masks and respirators have a medical-grade adhesive that conforms around the nose, therefore these strapless NIOSH-approved N95 masks do not contain any metals, metallic fibers, antimicrobial coatings, or metallic nanoparticles. Further, these strapless metal-free masks and respirators are safe for use in MRI procedures.

Please feel free to contact us with any additional questions.

Sincerely,

A handwritten signature in purple ink that reads "Kirstie Conklin".

Kirstie Conklin
Director of Marketing
Alliant Family of Companies
kconklin@allianthealthcare.com