



MRI (Magnetic Resonance Imaging) Safety Declaration

June 12, 2024

This is a Safety Declaration for anyone receiving a Magnetic Resonance Imaging exam (MRI). The U.S. Food and Drug Administration (FDA) is informing healthcare providers and patients that they may be injured if they wear face masks (surgical masks, nonsurgical masks, and respirators) with metal parts and/or metal coatings during a Magnetic Resonance Imaging (MRI) exam. The metal parts attached to many masks, including nose pieces or nose clips, wires, nanoparticles (ultrafine particles), or antimicrobial coating that may contain metal (such as silver or copper), can become hot, burning the patient during an MRI.

However, the FDA recommends that patients wear face masks without metal during MRIs. Alliant Biotech masks and respirators have a medical-grade adhesive that conforms around the nose. Therefore, these strapless NIOSH-approved N95 masks do not contain any metals, metallic fibers, antimicrobial coatings, or metallic nanoparticles. Further, these strapless metal-free masks and respirators are safe for use in MRI procedures.

Please feel free to contact us with any additional questions.

Sincerely,

A handwritten signature in black ink that reads "Kelli Jonas". The signature is written in a cursive, flowing style.

Kelli Jonas
Vice President, QA/RA
Alliant Healthcare Group
kjonas@medisurge.com