

This document outlines registered devices that can provide an adjunctive delivery option for surgeons desiring percutaneous targeted delivery and controlled access to bone defects, subchondral defects, lesions and cysts.

MIXING & DELIVERY SYSTEM

A delivery system that provides complete mixing of Mg OSTEOCRETE and excellent visibility for precise graft placement.

| Part Number | Description |
|---------------|--------------------------|
| 44-000-00-AUX | Mixing & Delivery System |



DELIVERY CANNULA

A delivery cannula that provides complete command over placement while an ergonomic T-handle allows for easy control. The universal luer lock connector ensures an easy fit for any size delivery syringe and the radiolucency of Mg OSTEOCRETE provides excellent visibility for precise graft placement.

| Part Number | Description |
|-------------|--|
| 74174-15M | Mg OSTEOCRETE™ Bone Marrow Aspiration Needle 8 Ga |
| 74066-01M | Mg OSTEOINJECT™ Bone Marrow Aspiration Needle 11 Ga |

8 Ga x 150mm, Open Tip, End Port Delivery 11 Ga x 110mm, Open Tip, End Port Delivery



* Distributed by Bone Solutions

DRILLABLE CANNULAS

Drillable cannulas with streamlined design for more versatile delivery options.

| Part Number | Description |
|-------------|---|
| 74358-01M | Mg OSTEONJECT™ 11 Ga Closed Tip Drillable Delivery Cannula |
| 74357-01M | Mg OSTEONJECT™ 11 Ga Open Tip Drillable Delivery Cannula |

* Distributed by Bone Solutions



SYRINGE PACK

Enables transport of prepared Mg OSTEONJECT from mixing syringe to smaller injectable quantity. Consists of five 1cc syringes and two connectors that attach to the mixing syringe.

| Part Number | Description |
|-------------|------------------------------------|
| K19-00059 | Mg OSTEONJECT™ Syringe Pack |

* Distributed by Bone Solutions



BEAD MAT KIT

Designed for targeted of Mg OSTEOCRETE to osseous defects and for dead space management.

| Part Number | Description |
|-------------|--------------------|
| RM44-050 | Bead Mat & Scraper |

